SleepImage 用户手册

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SleepImage® System 2.28 版

MyCardio LLC
3200 E Cherry Creek South Dr.
Suite 540
Denver, CO 80209, USA
support@sleepimage.com
www.sleepimage.com

SleepImage System 是一种医疗器械

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The screenshot shows the 'Account Settings' section of the SleepImage system. At the top, there is a navigation bar with links for Patients, Clinicians, Profile, Settings (which is highlighted in orange), and Help. Below the navigation bar, the title 'Account Settings' is displayed. A message states 'Account settings have been moved' and instructs users to manage their password and 2FA settings by clicking a button below. There is a 'Edit Account Settings' button. On the right, there is a 'Automatic Patient Sharing' section with a 'Email Address' input field and an 'Update' button. A horizontal line separates this from the 'Site-Wide Settings' section.

Site-Wide Settings

These settings will affect all your users

Allow Clinicians to view all patients.


When this feature is turned *on*, clinicians are able to view all the patients entered into your account and access their study data.

When this feature is *off*, you assign which clinicians are allowed to view which patients.

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
Enter Country or Time Zone name to set the Time Zone for your location

	America/Denver
--	----------------

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The screenshot shows the SleepImage Store interface. At the top, there's a dark header bar with the SleepImage logo, a navigation menu (Patients, Clinicians, Profile, Settings, Credits, Help), and user information (Role: Admin, SleepImage Test Account). Below the header is a "Patient List" section with a table showing patient details like name, gender, and status. A "Download all studies (CSV)" button is visible at the bottom right of this section.

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The screenshot shows the SleepImage My Account page. It includes sections for "CURRENT SUBSCRIPTION" (Business plan, \$11,776.86 per year, with a "View details" link and a "Cancel subscription" button), "PAYMENT METHOD" (Visa card ending in 4242, with a "link" and "..." button, and a "+ Add payment method" link), "BILLING INFORMATION" (Billing address: 80209 US, with a "Update information" link), and "INVOICE HISTORY" (a table showing a single entry for Mar 13, 2025, amount \$11,776.86, status Paid, and category Business). A "Return to SleepImage" link is also present.

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

The screenshot shows the 'Add New Clinician' dialog box. It includes fields for the clinician's name ('democlinician'), email ('democlinician@sleepimage.com'), and a checked checkbox for sending an email to set a password. Below these are fields for first name ('Demo'), last name ('Clinician'), middle initial ('optional'), specialty ('Family Medicine'), and language ('English (U.S.)'). A dropdown for time zone ('Enter Country or Time Zone name to set the Time Zone for Clinician location') shows 'America/Denver'. At the bottom are 'Add' and 'Cancel' buttons.

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The screenshot shows the 'Edit Clinician' dialog box. It includes fields for the clinician's name ('democlinician'), email ('democlinician@sleepimage.com'), and password ('New Password' and 'Confirm New Password'). Below these are fields for first name ('Demo'), last name ('Clinician'), middle initial ('optional'), specialty ('Family Medicine'), and language ('English (U.S.)'). A dropdown for time zone ('Enter Country or Time Zone name to set the Time Zone for Clinician location') shows 'America/Denver'. At the bottom are 'Deactivate User', 'Save', and 'Cancel' buttons.

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The screenshot shows the 'Add New Patient' window with the following fields filled:

- Demographic:**
 - Demo:** Demo
 - Patient:** Patient
 - Middle Initial (optional):** (empty)
 - Patient Number (optional):** 01/01/1970
 - Gender:** Male (radio button selected)
 - Email (optional):** (empty)
 - Language:** English (U.S.)
- Time Zone:** America/Denver
- Current Study Count:** 3
- Enable study limit:** checked
- Buttons:** Add (dark blue), Cancel (orange)

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Current Study Count:

Patients are allowed unlimited studies when this setting is off. Only studies ≥ 4 hours will count against the limit. After this limit is reached, additional studies will not be automatically processed. Studies can be processed beyond this limit with manual approval.

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Edit Patient

Demographic Vitals Sleep Complaints Medical History Medications & Treatments

Demo

Patient

Middle Initial (optional)

Patient Number (optional)

01/01/1970

Male Female Other

Email (optional)

English (U.S.)

Enter Country or Time Zone name to set the Time Zone for Patient location

America/Denver

Current Study Count:
Patients are allowed unlimited studies when this setting is off. Only studies ≥ 4 hours will count against the limit. After this limit is reached, additional studies will not be automatically processed. Studies can be processed beyond this limit with manual approval.

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

The screenshot shows the 'Account Settings' section of the SleepImage web interface. On the left, under 'Site-Wide Settings', there is a box containing the 'Allow Clinicians to view all patients.' checkbox. This checkbox is currently unchecked. To its right, under 'Authorized Connection', is a section for 'SleepImage API' which also includes an unchecked checkbox. At the bottom of the page, there is an 'Automatic Patient Sharing' section with a 'Email Address' input field and an 'Update' button.

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The screenshot shows the SleepImage web interface. On the left, the 'Patient List' table displays three patients: Green, Michael; Johnson, Emily; and Smith, John. The row for Johnson, Emily is highlighted. A modal dialog box titled 'Share Patient' is open over the patient list. It contains a recipient email input field with 'email@email.com' and a 'Close' button at the bottom. The background shows the SleepImage header with 'Sign Out' and 'Role: Clinician | 2.23 Testing Clinic'.

Name	Last Study	Studies >4hrs / T...	Study Limit	Patient Num...
Green, Michael	Pending	1 / 1	1	
Johnson, Emily	08/10/2024	2 / 3	3	
Smith, John	04/06/2024	3 / 3		

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The screenshot shows the SleepImage portal interface. At the top, there's a navigation bar with links for Patients, Profile, Settings, Help, and Sign Out. Below this, on the left, is the 'Account Settings' sidebar with sections for Change Password, 2-Factor Authentication (2FA), and a general settings section. The main content area is titled 'Add template'. It includes fields for Date (Sep 18, 2024), Template name, Template type, Status, Private note, and Template content. There are 'Save' and 'Cancel' buttons at the bottom. In the background, partially visible, is a table of interpretation templates with columns for Title, Type, Date Added, Notes, and Content.

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The screenshot shows two dialog boxes. The left one is titled 'Edit template' and contains fields for Date (Sep 3, 2024), Template name (Admin - Example Impression), Template type (Impression), Status (Active), and a large 'Template content' text area. The right dialog box is titled 'Confirm Remove Template' and asks 'Are you sure you want to remove this template?' with 'Yes' and 'No' buttons. Above the 'Yes' button, it shows the template details: Admin - Example Impression (09/03/2024).

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Upload Studies

Last upload Sep 17, 2024 9:48 AM was successful.

V223 Testing, Verification (01/01/1970)
Upload By: V2.23 Testing

Processing is complete.

Select a checkbox next to each study to either Accept or Reject the study for report generation.

- Each Accepted study will generate a Report to the Patient's Profile.
- Each Rejected study will not generate a Report and is not added to the Patient's Profile.

Once a selection is made for each study, click Confirm to continue, or Cancel to delay selection and close the window.

Accept	Reject	Start Date Time	Duration (HH:MM)	Status
<input type="checkbox"/>	<input type="checkbox"/>	04/04/2024 01:01:34 AM	11:42	Success

1 studies have been uploaded.

- 0 are being accepted for report generation.
- 0 are being rejected.

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SleepImage

Patients Clinicians Profile Settings Credits Help

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
Sleep Complaints:

Select all that apply

Excessive Daytime Sleepiness Irregular / Paused Breathing During Sleep
 Snoring Bruxism / Teeth Grinding
 Morning Headaches Waking Up Gasping Or Choking
 Difficulty Falling Or Maintaining Sleep

This is the first admission for this 56 year old woman, Convey the acute or chronic nature of the problem and who states she was in her usual state of good health until establish a chronology. one week prior to admission. At that time she noticed the abrupt onset (over a few seconds to a minute) of chest pain onset

Character Count: 961 / 5000

Medications Used:

1. Pepcid 2. Tylenol [Add More](#)

Treatment Device Used:

1. Helmet [Add More](#)

Diagnosis:

Select all that apply

G47.33 Obstructive Sleep Apnea G47.31 Primary Central Sleep Apnea
 G47.32 High Altitude Periodic Breathing G47.9 Sleep Disorder, Unspecified
 G47.10 Hypersomnia, Unspecified G47.00 Insomnia Disorder
 R06.3 Periodic Breathing R06.83 Snoring
 R53.83 Other Fatigue Other

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If you would like to send email notification(s), please enter the email addresses below, separated by a comma. If this note is for a shared patient for which you are the recipient, the patient owner's (the sender) email address will be automatically populated in the field

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Medical Data Changed

You have unsaved Patient Medical Data.

Would you like to save these changes to the server?

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The screenshot shows the 'Study Details' section of a software interface. At the top right, it says 'Role: Admin | Demo Clinic'. Below that is a 'Selected Page' dropdown set to 'Reports'. To the right of the dropdown is a 'Raw Data' button with a blue circular icon. At the bottom left is a 'Sleep Report' summary box containing patient information: 'Patient Name: Demo Patient', 'Patient Id: YSN2PP', and 'Date Of Birth: 01/01/1970 (54 years)'. The main area is titled 'Sleep Report' and shows a timeline from 'Apr 3, 2024 11:48 PM - 7:02 AM' to 'Demo Clinic'. It includes sections for 'CFC Freq (Hz)', 'Stable NREM', 'Unstable NREM', 'REM/WAKE', 'Fragmentation', 'Periodicity', 'CVHR', 'Desat.', 'AHI', 'Events Cent.', 'SpO2 (%)', and 'Actigraphy Signal'. A vertical orange arrow points to the '10:35 PM' mark on the timeline.

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版权声明 ©