SleepImage Ring

Model: PO4

User Manual

This device uses Bluetooth to communicate with the SleepImage Mobile App. The device must be paired with the App to collect and transfer data collected during sleep to the SleepImage System, Software as a Medical Device.

Download App

Name: [SleepImage](#)

iOS: App Store

Android: Google Play

Please ensure App is updated to the latest version prior to each use.

1. Introduction

1.1 Intended use

This Pulse Oximeter is intended to be used for measuring, displaying and storing adult's pulse oxygen saturation (SpO2), pulse rate of adults in home or healthcare facilities environment. (EU countries)

This product is intended to be used with the SleepImage System for recording plethysmograph, oxygen saturation (SpO2) and actigraphy data that is transferred via Bluetooth to the SleepImage App during recording sessions. (USA)

The Pulse Oximeter is a pulse oximeter indicated for use in measuring, displaying, storing and transmitting functional oxygen saturation of arterial hemoglobin (SpO2) and pulse rate for adult patients. It is intended for spot-check and/ or continuous data storage, and not continuous monitoring. It can be used in sleep labs, long-term care, hospitals (except the operating room, ICU) and home use. (Australia & Canada)

Note: Data storage is not developed in the current version.

1.2 Warnings and Cautions

- DO NOT squeeze the sensor or apply excessive force.

- Do not use the Ring during an MRI examination.
- Do not use the Ring with a defibrillator.
- Do not store the Ring in the following locations: direct sunlight;


high temperatures, high levels of moisture, or heavy contamination; locations near to sources of water or fire; or locations subject to strong electromagnetic influences.

- Do not use the Ring in a combustible environment.
- Never submerge the Ring in water or other liquids.
- Do not clean Ring with acetone or other volatile solutions.
- Do not drop the Ring or subject it to strong impact.
- The Ring and accessories are non-sterile.
- Do not place Ring in pressure vessels or gas sterilization devices.
- Do not dismantle the Ring.
- Consult your doctor immediately if you experience symptoms that could indicate acute disease.
- Do not self-diagnose or self-medicate solely on the basis of data displayed by the Ring.
- Only use the charging cable specified in this manual.
- Wearing the Ring during a normal sleep duration is generally safe but may cause skin irritation.
- Do not open the device cover without authorization. The cover should only be opened by a qualified service personnel.
- The device shall only be maintained by qualified professionals.
- The manufacturer shall provide the service personnel with circuit diagrams, component part lists, descriptions, calibration instructions, or other information that will assist service personnel to repair the device.

1.3 Included in Box

Device, User Manual, Warranty Card, Charging Cable.

2. Display Overview

3. Using the Ring and App

3.1 Charging

Charge battery before using by connecting Ring to computer USB or USB charging adapter with the included USB cable. When fully charged, Ring will power off automatically.

3.2 Display

Touch the SleepImage logo on the display to switch between displaying (1) time & battery level, (2) serial number and (3) Pulse/SpO2 and unavailable symbol.

3.3 Screen Wake up

The screen turns off automatically during recordings; you can

touch the SleepImage logo on the display to wake up the screen.

3.4 POWER ON/OFF

Device will turn on automatically by placing a finger in the sensor and turn off automatically by removing a finger from the sensor.

3.5 SleepImage App

The Ring uses Bluetooth to communicate with the SleepImage Mobile App. The Ring must be paired with the App to collect and transfer data to SleepImage.

Note: DO NOT PAIR in the settings of your phone. Follow the Mobile App Instructions for Use to download the App and pair the Ring before attempting data collection.

3.6 Start recording

Wear Ring on thumb or index finger. Ring will turn on automatically. **Start the recording in the App to begin data collection.**

Ring fit and other notices:

- Wear Ring on thumb or index finger.
- A proper fit is essential.** If the Ring is too tight for thumb and index fingers, try an alternate finger, but avoid the middle finger. A loose fit may lead to poor data collection.
- Avoid excessive motion during recording.
- Avoid strong ambient light during recording.
- If the Ring is removed briefly during a recording session, simply put the Ring back on to resume the recording.

3.7 Stop recording

After waking, Stop the recording in the App and remove the Ring which will automatically turn itself off.

3.8 Data Upload

After the recording is stopped in the App, the recording will upload to SleepImage. If study does not upload, check Troubleshooting in the SleepImage App instructions on www.sleepimage.com.

3.9 Unavailable Symbol

When this symbol displays on device screen, it indicates the readings is unavailable, which may be caused by:

- Excessive movement;**
- Poor signal;**
- Finger is too cold.**

4. Maintenance

4.1 Time & Date

Device time automatically syncs from your mobile device.

4.2 Cleaning

Use a soft cloth lightly moistened with water or isopropyl alcohol to clean the device surface.

4.3 Recycling

Disposal of the device and its accessories must be done in accordance with local regulations for recycling.

5. Troubleshooting

Problem	Possible Cause	Possible Solution
Device does not turn on or no response	Battery may be low.	Charge battery and try again.
	Device may be damaged.	Please contact your local distributor.
The app cannot find or pair with the device	The device is off.	Turn on device.
	Bluetooth cannot work without location permission (Android)	Allow location access.
Only one Light Emitter on the ring turns red.	The Bluetooth of your phone is off.	Turn on the Bluetooth in the phone.
	This is normal, the O2Ring only has one light emitter.	No need to worry about it.
The device screen displays "Error 1"	Errors occur during data analysis.	Connect the power supply and plug in the charging cable, press and hold the touch key for 3s to reset the hardware.
	Your mobile device may not meet the minimum specifications required. (See www.sleepimage.com for information).	

6. Guide to Symbols

Symbol	Description
	Type BF-Applied Part
	Manufacturer
	Date of manufacture
	Follow Instructions for Use.
	No alarm system
	MRI unsafe. Presents hazards in all MR environments as device contains strongly ferromagnetic materials.
	Conforms with Federal Communications Commission standards for electromagnetic safety.
	Declaration of compliance with EU Medical Device Directive.
	Authorized representative in the European community.
	UKCA marking.
	Authorized Representative in the United Kingdom.
	Non-ionizing radiation.
	This device contains license-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada’s license-exempt RSS(s).
	Against ingress of solid foreign objects ≥ 12.5mm diameter, splashing.
	Medical device
	Serial number
	Temperature limitation
	Humidity limitation
	Atmospheric pressure limitation
	Indicate separate collection for electrical and electronic equipment (WEEE).
	This product complies with verpackG.
	Our products and packaging can be recycled, don’t throw them away! Find where to drop them off on the www.quefairedemesdechets.fr site (Only applicable for French market)

7. Specifications

Environmental	Operating	Storage
Temperature	5 to 40°C	-25 to 70°C
Relative humidity (noncondensing)	10% to 95%	10% to 95%
Barometric	700 to 1060hPa	700 to 1060hPa
Protection against electric shock	Internally powered equipment	
Degree protection against electrical shock	Type BF	
Electro-magnetic compatibility	Group I, Class B	
Degree of dust & water resistance	IP22	
Weight	15 g	
Size	38mm×30mm×38 mm	
Battery	3.7Vd.c., Rechargeable Lithium-polymer	
Charge time	2-3 hours	
Battery life	14 hours for typical use	
Wireless	Bluetooth 4.0 BLE	
Oxygen level range	70% to 100%	
SpO2 Accuracy (Arms)	80-100%:±2%, 70-79%:±3%	
Pulse Rate range	30 to 250 bpm	
Pulse Rate accuracy	±2 bpm or ±2%, whichever is greater	
Recorded parameters	Oxygen level, Pulse Rate, Plethysmograph, Actigraphy	
Data storage	No data is stored on the Ring	
Frequency range	2.402 – 2.480 GHz	
Max RF power	-10 dBm	
Expected service life	3 Years	
Wavelength/Max emission power	660nm/940nm, 0.8mW/1.2mW	

8. FCC Statement

FCC Warning:
FCC ID: 2ADXX-8326
Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and
(2) this device must accept any interference received, including interference that may cause undesired operation.
Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
-Reorient or relocate the receiving antenna.
-Increase the separation between the equipment and receiver.
-Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
-Consult the dealer or an experienced radio/TV technician for help.
The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

9. IC Caution

IC:29845-S001
This device contains license-exempt transmitter(s)/ receiver(s) that comply with Innovation, Science and Economic Development Canada’s license-exempt RSS(s). Operation is subject to the following two conditions:
(1) This device may not cause interference.
(2) This device must accept any interference, including interference that may cause undesired operation of the device.
Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

10. Electromagnetic Compatibility

The device meets the requirements of IEC 60601-1-2.

Shenzhen Viatom Technology Co., Ltd.
901, Building West, Lepu Tower, No.66 Xingke Road,
Xili Community, Xili Street, Nanshan District,
Shenzhen, 518055, Guangdong, P.R. China
www.viatomtech.com

EC

REP

MedNet EC-REP GmbH
Borkstrasse 10, 48163 Muenster, Germany
Tel: +49 251 32266-0 Fax: +49 251 32266-22
Email: contact@mednet-ecrep.com

UK

REP

MediMap Ltd
2 The Drift, Thurston, Suffolk IP31 3RT, United Kingdom
Tel: +49 251 32266-0 Fax: +49 251 32266-22
Email: contact@mednet-ecrep.com

Australia Sponsor: SHARE INFO PTY LTD
Add: 4 Allnutt ct, Cheltenham, Melbourne, VIC 3192,
Australia

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Trade Name: SleepImage Ring Device
Product Name: Pulse Oximeter
Version: V1.0 Date: February 6, 2026 Model: PO4 PN: 255-03413-01
Distributed by MyCardio LLC, 3200 E Cherry Creek South Dr., Suite 540, Denver, CO 80209, USA
Rev 6 Issued February 6, 2026
SleepImage Website: sleepimage.com
Contact SleepImage: support@sleepimage.com
Manufacturer Website: www.viatomtech.com
Contact Manufacturer: service@viatomtech.com

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