SleepImage Ring

Model: PO2

User Manual

The SleepImage Ring (Ring) is a version of the Viatom model PO2 and uses Bluetooth to communicate with the SleepImage Mobile App (App). The Ring must be paired with the SleepImage App to collect and transfer data to the Sleepimage System, a Software as a Medical Device (SaMD).

Download the app by searching: <https://sleepimage.com/patient-resources/>

The Ring displays Pulse Rate and SpO₂ data while in use. It should not be used to directly diagnose or treat any medical condition, and it does not store any data.

1. Introduction

1.1. Intended use

The Pulse Oximeter is a non- invasive device indicated for use in measuring, displaying, storing and transmitting functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate for adult patients. It is intended for spot-check and/or continuous data collection, and not continuous monitoring. It can be used in sleep labs, long-term care, hospitals and home use. (FDA, CA)

This Pulse Oximeter is intended to be used for measuring, displaying and storing of pulse oxygen saturation (SpO₂), pulse rate of adults in home or healthcare facilities environment. (CE, AU)

Intended users: adults.

1.2. Name and Model

Name: Pulse Oximeter

Model:PO2

1.3. Warnings and Cautions

- Federal law restricts this device to sale by or on the order of a physician.
- DO NOT squeeze the sensor part or apply excessive force on it.

- DO NOT use this device during MRI examination.
- DO NOT use this device with a defibrillator.
- DO NOT use this device if you have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device. Such use could cause electrical shock, burns, electrical interference, or death.
- DO NOT store the device in the following locations: locations exposed to direct sunlight, lint, dust, high temperatures, high humidity, heavy contamination, locations near sources of water or fire, or locations subject to strong electromagnetic influences.
- DO NOT use the device in a combustible environment, including oxygen-enriched environments.
- NEVER submerge the device in water or other liquids.
- The device and its accessories are provided non-sterile.
- DO NOT clean the device with acetone or other volatile solutions.
- DO NOT drop the device or subject it to strong impact.
- DO NOT place the device in pressure vessels or gas sterilization devices.
- DO NOT dismantle, modify, or open the device without authorization, as this may cause damage, malfunction, or affect device performance. The device cover should only be opened by qualified service personnel.
- The device shall only be maintained by qualified professionals.
- The manufacturer shall provide service personnel with circuit diagrams, component part lists, descriptions, calibration instructions, and other information necessary for device repair.
- Consult your doctor immediately if you experience symptoms that may indicate acute disease.
- DO NOT self-diagnose or self-medicate based solely on the data displayed by this device. Do not start taking new medication or change the type and/or dosage of any existing medication without prior approval from your doctor.
- Use only cables, sensors, and other accessories specified in this manual.
- Prolonged continuous recording or monitoring may increase the risk of undesirable skin changes, such as irritation, redness, blistering, or burns. Check the wearing position or sensor application site every 6-8 hours.
- The biocompatibility testing has been performed on the materials in contact with the user in accordance with ISO 10993.
- Do not wear the device on fingers with edema, fragile tissue, or damaged skin.

- Check the SpO₂ sensor and cable before use. Do not use a damaged SpO₂ sensor.
- Check the SpO₂ sensor application site every 6-8 hours to verify sensor positioning, circulation, and skin sensitivity. Patient sensitivity may vary depending on medical condition or skin condition. For patients with poor peripheral blood circulation or sensitive skin, inspect the sensor site more frequently.
- The functional tester cannot be used to assess the accuracy of the SpO₂ sensor or the device.
- The device has no alarm system.
- Follow local laws and regulations when disposing of the device and accessories.
- Do not maintain the device while it is charging.
- Keep the cable away from children. The cable may cause strangulation.
- Keep the device out of reach of pets, pests, and children.
- The PULSE OXIMETER EQUIPMENT is calibrated to display FUNCTIONAL OXYGEN SATURATION.
- This device is designed to determine the arterial oxygen saturation percentage of functional hemoglobin. Factors that may affect pulse oximeter performance or measurement accuracy include, but are not limited to:
 - Excessive ambient light;
 - Excessive motion;
 - Electrosurgical interference;
 - Blood flow restrictors (such as arterial catheters, blood pressure cuffs, infusion lines, etc.);
 - Moisture in the sensor;
 - Improper sensor application;
 - Incorrect sensor type;
 - Poor pulse quality;
 - Venous pulsations;
 - Anemia or low hemoglobin concentrations;
 - Cardiogreen and other intravascular dyes;
 - Carboxyhemoglobin;
 - Methemoglobin;
 - Dysfunctional hemoglobin.

- DO NOT use this device during low perfusion conditions.
- DO NOT use this device for SpO₂ measurement while moving and/or under low perfusion conditions, as this may affect the accuracy of the SpO₂ measurement.
- When the signal is incomplete (for example, due to excessive signal noise, poor signal quality, or signal loss), the blood oxygen and pulse rate values may become invalid. The displayed values will disappear and the display will show "--".
- The device must be allowed to stabilize at 20 °C for 1 hour after removal from the minimum or maximum storage temperature environment before normal use.
- Although the device complies with electromagnetic disturbance immunity requirements, under certain conditions (such as strong electromagnetic interference from diathermy, electrosurgery, radiogram equipment, television, radio communication equipment, RFID systems, or security systems), device performance cannot be guaranteed. If abnormal indicators are displayed (such as pulse strength, pulse symbol, or pulse plethysmogram), the device may be affected by electromagnetic disturbance and measurement results may be inaccurate.
- DO NOT operate this device in environments with strong electromagnetic interference, such as MRI scanning rooms or areas containing high-frequency electrosurgical equipment.
- Only use power adapters specified in this manual. Power adapters should be 5 V, 1 A or higher and comply with EN/IEC 60601-1 requirements.
- Laptops or external equipment connected to the device should comply with IEC 60950-1 and/or IEC 62368-1 requirements.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

1.4. Clinical Benefit

Clinical Benefits: Providing the measurement of both the pulse oxygen saturation (SpO₂) and pulse rate data for adult patients in homes and medical clinics.

1.5. Contraindications

No contraindications.

1.6. Guide to Symbols

Symbol	Description
	Manufacturer

Symbol	Description
	Date of manufacture
	Serial number
	Unique device identifier
	Indicates that the product should not be discarded as unsorted waste but must be sent to separate collection facilities for recovery and recycling.
	Refer to instruction manual
	Type BF Applied Part
	No alarm system
	MR unsafe
	Indicates that the product is protected against solid foreign objects of 12,5 mm Ø and greater; and protected against vertically falling water drops when enclosure tilted up to 15°.
	Caution
	General warning sign
	Indicates that the product complies with the EU Medical Device Regulations(Regulation (EU) 2017/745)
	Conforms with Federal Communications Commission standards for electromagnetic safety.
	This device contains license-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's license-exempt RSS(s).
	Medical device
	Authorized representative in the European community
	UKCA marking
	UK Responsible Person
	Non-ionizing radiation
	Indicates that the marked item or its material is part of a recovery or recycling process.
	Our products and packaging can be recycled, don't throw them away! Find where to drop them off on the www.quefairedemesdechets.fr site (Only applicable for French market).
	Temperature limit
	Humidity limitation
	Atmospheric pressure limitation

1.7. Included in Box

SleepImage Ring, Quick Start Guide, Warranty Card, Charging Cable.

2. Using the Ring and App

2.1 Charging

Charge battery before using by connecting Ring to computer USB or USB charging adapter with the included USB cable. A fully charge will need 2-3 hours. When fully charged, Ring will power off automatically.

Note: The device cannot be used during charging, and if choosing a third-party charging adaptor (Class II), select one that complies with IEC60950 or IEC60601-1.

2.2 Download the App

SleepImage Mobile App

App Name: SleepImage

iOS: App Store

Android: Google Play

NOTE:

- Ensure that the App is updated to the latest version prior to each use.


2.3 Pair the Ring with the App

The Ring uses Bluetooth to communicate with the SleepImage Mobile App. The Ring MUST be paired with the App to collect and transfer data to Sleepimage.

NOTE:

- Only pair the Ring to your phone in the App. DO NOT pair in phone settings and follow the Mobile App Instructions for Use.

2.4 Ring Display

Touch the SleepImage logo on the display to switch between displaying (1) time & battery level and (2) Pulse/SpO₂ readings or unavailable symbol.

2.5 Unavailable Symbol

This symbol on the Ring display indicates that readings are temporarily unavailable, which may be caused by:

- Excessive movement;
- Poor signal;
- Finger being too cold;
- Ring initiating at the start of a recording.
- Low Perfusion


2.6 Screen Wake up

The screen turns off automatically during recordings. Touch SleepImage logo on the display to wake up screen.

2.7 POWER ON/OFF

The Ring will turn on automatically when worn and turn off when removed from finger.

2.8 Start recording

- Wear the device on thumb finger. Move the device along the finger.
- Device will turn on automatically. After a few seconds, the device will begin to work.

NOTE:

- Proper fit is essential. If the Ring is too tight for thumb or index finger, try an alternate finger, but avoid the middle finger.
- Loose fit may lead to poor data collection.
- Avoid excessive motion during recording.
- Avoid strong ambient light during recording.
- If the Ring is removed briefly during a recording session, simply put the Ring back on to resume the recording.
- DO NOT stop recording in the App during the night if you wake up and temporarily take the Ring off.

SpO2 measurement principle:

The Pulse Oximeter is a lightweight, portable health oximeter for use in the home or in healthcare facilities. SpO2 measurement technology is based on developed photoelectron method, the circuit design and calculation software were developed by Shenzhen Viatom Technology Co., Ltd. The SpO2 sensor receives the optical signal from the red light and infra-red light through the finger. Insert the finger into the oximeter, there are two emitting tubes (red light diodes and infrared diodes) located on the inner upside of the sensor and they can emit red light and infrared; There is the receiving end located on the inner downside of the sensor, and it can transmit the red light and infrared into the pulse signal through finger. The MCU receives the pulse signal, gets the frequency signal by counting, processes its digital signal, and finally gets the measured SpO2 value. The PR is averagely calculated by above peak intervals of PR waveform.

2.9 Stop Recording

At the end of your sleep test, stop the recording in the App and remove the Ring from your finger which will automatically turn itself off.

2.10 Data Upload

After the recording is stopped in the App, the recording will upload to

SleepImage. If the study does not upload, check Troubleshooting in the Sleepimage App instructions, or find further information on www.sleepimage.com.

2.11 Cybersecurity

- The device can transmit clinical data (SpO₂ data) from the subject device to a wireless receiver for storage and display. The common receiver is APP software. So it is needed to install international mainstream anti-virus software such as AVG, Avast, McAfee, etc. In the running environment, even update the virus upgrade package in time and scan for viruses. For IOS system, the APP can run without antivirus software and firewalls because of the closed system and access control of the system. The end user needs to select the device ID before they want to conduct Bluetooth connection.

- In order to support the cybersecurity of this device, we provide information about the risk management with regards to cybersecurity which includes the assets, vulnerabilities, threats and controls. If you have any problems, please contact us.

Identification	Risk Description	Impact
Threats	CR1: Malicious attack such as Counterfeit command	Maybe the device failed to use and the treatment is delayed.
	CR2: Malicious attack such as Violent login attempt / password guessing attack	Maybe the treatment is delayed.
	CR3: Sensitive information including health data and the user information leakage	Maybe Tampering and Embezzlement of the user' information.
	CR4: Wrong operation of critical and/or sensitive information by authorized users	Maybe Tampering and Embezzlement of the user' information.
	CR5: Improper Authority Setting	Maybe the device failed to use.
Vulnerabilities	CR6: Access without authentication	Maybe the device failed to use.
	CR7: Sensitive information with encryption	Maybe the user' data is lost or stolen.
	CR8: Sensitive information with integrity	Maybe the user' data is lost or stolen.
Assets	CR9: Information Assets: Patient data; Device settings; Health data; Connection ID	Maybe the health data is revised or lost.
	CR10: Physical Assets :User interface; Network interface; Device assets (Operating system, Software etc.)	Maybe the device failed to use.

- The device only includes the embedded software. The firmware is directly coded into the integrated system with on circuit unchangeable package. There is no software transferring during production or distribution. The manufacturer does the burning and upgrading of the subject device software and there is no case of being attacked by malware. Then installing antivirus software and firewalls in the device is unnecessary.

- SleepImage regularly checks whether the released app software is the latest version, whether it can be downloaded and installed, and whether the software is attacked (MD5 hash value verification) to ensure software integrity and authenticity.

- All updates and patches will be released to customer through distributor. A release note or letter will be issued and the version number will be changed accordingly.

- Device software upgrade can be done remotely through our APP. When the device is connected to the APP, the APP will check the software version of the device, and when a new available version is found, an upgrade button will pop up and the user will authorise the upgrade operation.

- Before discarding the device, it needs to be ensured that all health

measurements and user account data stored in the device and the app are secure and confidential (e.g. deleted).

3. Maintenance

3.1 Time & Date


After connection with App, device time will sync from your phone time automatically.

3.2 Cleaning

Use a soft cloth moistened with water or alcohol to clean the device surface.

3.3 Disinfection

Use 75% alcohol wipes to clean the screen of the device and the crevices and around the edges of the device, and then let it air dry. As shown in the picture below:

Use 75% alcohol wipes to clean the SpO₂ sensor carefully, and then let it air dry. As shown in the picture below:

CAUTION: To avoid permanent damage to the device, do not use petrol, thinners or similar solvent.

CAUTION: After cleaning, check the surface and, the interiors of the Finger Nail Position and especially the Lighting Position under using a magnifying glass or sufficient light.

CAUTION: If dirt is found, please repeat the above steps using a fresh wipe.

CAUTION: If you find that the surface or inside of the finger cover is cracked or the product cannot be turned on, please contact us immediately.

The device is a non-sterile medical device and does not contain any sterile or degradable component thus the device is not subject to the shelf life requirements.

In order to ensure its long service life, please pay attention to the maintenance.

4. Troubleshooting

Problem	Possible Cause	Possible Solution
Device does not turn on or no response	Battery may be low.	Charge battery and try again.
	Device might be damaged.	Please contact your local distributor.
	Software exception	Keep device in charging, touch the key for 8 seconds.
The app cannot find the device	The Bluetooth of your phone is off.	Turn on the Bluetooth in the phone.
	The device Bluetooth is off.	Turn on device
	For Android, Bluetooth cannot work without location permission	Allow location access
Only one Light Emitter on the ring turns red.	This is normal, the device only has one light emitter.	No need to worry about it.
The device screen displays “Error 1”.	Errors occur during data analysis.	Connect the power supply and plug in the charging cable, press and hold the touch key for 3s to reset the hardware.
Your mobile device may not meet the minimum specifications required. (See www.sleepimage.com for information).		

5. Specifications

Environmental	Operating	Storage
Temperature	5 to 40°C	-10 to 45°C
Relative humidity (non-condensing)	10% to 80%	10% to 93%
Barometric	700 to 1060hPa	700 to 1060hPa
Protection against electric shock	Internally powered equipment	
Degree protection against electrical shock	Type BF	
Electro-magnetic compatibility	Group I, Class B	
Degree of dust & water resistance	IP22	
Weight	15 g	
Size	38mm×30mm×38 mm	
Battery	3.7Vd.c., Rechargeable Lithium-polymer	
Charge time	2-3 hours	
Charge requirement	DC 5V 1A	
Battery life	12-14 hours for typical use	
Wireless	Bluetooth 4.2 BLE	
Oxygen level display range	0% to 100%	
SpO ₂ Accuracy (Arms)	70-100%:±2%(Arms:1.77); 70-80%: ±3%; 80%-90%:±2%; 90%-100%:±2% 0%-69%: not defined	
Pulse Rate range	25 to 250 bpm	
Pulse Rate accuracy	±2 bpm or ±2%, whichever is greater	
A functional tester or SpO ₂ simulator (Model: PRONK TECHNOLOGIES OX-1) can be used to determine the pulse rate accuracy.		
Recorded parameters	Oxygen level, Pulse Rate	
Data storage	No data is stored on the Ring	
Data update period	<10s	
Frequency range	2.402-2.480GHz GFSK Modulation Adaptive Frequency Hopping (AFH)	
Wireless Quality of Service (QoS)	Transmission Distance: 1.5m Transmission Time: ≤10s Data integrity: 100%	
Network topology	Point-to-Point	
Band width	1Mbps	
Max RF power	-10 dBm	
Expected service life	3 years	
Mobile App for iOS	iOS 9.0 or above, iPhone 4s/ iPad 3 or above	
Mobile App for android	Android 5.0 or above, with <i>Bluetooth</i> 4.0 BLE	
Wavelength/Max emission power	660nm±5nm/940nm±10nm, 0.8mW/1.2mW	

6. FCC Statement

FCC Warning:

FCC ID: 2AD XK-8326

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed

to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

7. IC Caution

IC:29845-S001

This device contains license-exempt transmitter(s)/ receiver(s) that comply with Innovation, Science and Economic Development Canada's license-exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

8. Appendix EMC

The device meets the requirements of IEC60601-1-2:2020.

⚠ Warnings and Cautions

- Using accessories other than those specified in this manual may result in increased electromagnetic emission or decreased electromagnetic immunity of the equipment.
- The device or its components should not be used adjacent to or stacked with other equipment.
- The device needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.
- Other devices may interfere with this device even though they meet the requirements of CISPR.
- When the inputted signal is below the minimum amplitude provided in technical specifications, erroneous measurements could result.
- Portable and mobile communication equipment may affect the performance of this device.
- Other devices that have RF transmitter or source may affect this device (e.g. cell phones, PDAs, and PCs with wireless function).
- Electromagnetic fields avoid using the product near strong radio frequency (RF) signals or portable and/or mobile RF devices and/or specific RF emitters that are known sources of electromagnetic disturbance such as diathermy, electrocautery, RFID, security systems (e.g., electromagnetic anti-theft systems, and metal detectors). Interference from hidden RF emitters like RFID might cause packet loss and this will be visible as a “Poor Bluetooth Signal” message on the mobile application. Move away from the hidden RF emitter if this happens.

Essential Performance		
Essential Performance		The description of what the operator of the device can expect if the Essential Performance is lost or degraded due to electromagnetic disturbances
Oxygen level display range	0% to 100%	Please stop using the device immediately and contact the device manufacturer or distributor for service as soon as possible
SpO ₂ Accuracy	70%-100%: ±2% (Arms:1.77%)	
	70%-80%: ±3%	
	80%-90%: ±2%	
	90%-100%: ±2% 0%-69%: not defined	
Pulse Rate range	30 to 250 bpm	
Pulse Rate accuracy	±2 bpm or ±2%, whichever is greater	

Table 1

Guidance and manufacturer's declaration-electromagnetic emission	
The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Pulse Oximeter should assure that it is used in such	

an environment.		
Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The Pulse Oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Pulse Oximeter suitable for use in all establishments, including domestic establishments and those directly network that supplies buildings used for domestic purposes.
Harmonic emissions IEC61000-3-2	N/A	
Voltage fluctuations/flicker emissions IEC61000-3-3	N/A	

Table 2

Guidance and manufacturer's declaration-electromagnetic emission			
The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Pulse Oximeter should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment -guidance
Electrostatic discharge (ESD) IEC61000-4-2	±8 kV contact ±15kV air	±8 kV contact ±15kV air	Floors should be wood, concrete or ceramic tile. if floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/ burst IEC61000-4-4	±2kV for power Supply lines ±1 kV for input/output lines	N/A	N/A
Surge IEC 61000-4-5	±1kV line (s) to line(s) ±2kV line(s) to earth	N/A	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5% UT >95% dip in UT) for 0.5 cycle <40% UT (60% dip in UT) for 5 cycles <70% UT (30% dip in UT) for 25 cycles <5% UT >95% dip in UT) for 5 s	N/A	N/A
Power frequency (50Hz/60Hz) magnetic field IEC61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

Table 3

Guidance and manufacturer's declaration - electromagnetic immunity			
The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of The Pulse Oximeter should assure that it is used in such an electromagnetic environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of The Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance d=1.2√P d=1.2√P 80MHz to 800MHz d=2.3√P 800MHz to 2.5GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). b Field strengths from fixed RF transmitters, as determined by an
Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

			electromagnetic site survey ,a should be less than the compliance level in each frequency range .b Interference may occur in the vicinity of equipment marked with the following symbol.
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a: Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, and electromagnetic site survey should be considered. If the measured field strength in the location in which The Pulse Oximeter is used exceeds the applicable RF compliance level above, The Pulse Oximeter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating The Pulse Oximeter. b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.			

Table 4

Recommended separation distances between portable and mobile RF communication the equipment			
The Pulse Oximeter is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of The Pulse Oximeter can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Pulse Oximeter as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W(Watts)	Separation distance according to frequency of transmitter M(Meters)		
	150kHz to 80MHz d=1.2√P	80MHz to 800MHz d=1.2√P	80MHz to 2,5GHz d=2.3√P
0,01	N/A	0.12	0.23
0,1	N/A	0.38	0.73
1	N/A	1.2	2.3
10	N/A	3.8	7.3
100	N/A	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Copyright

This manual is written by our company and all rights reserved. Without our company's prior written consent, no part of this manual may be reproduced or copied in any form or method.

Illustration

All illustrations provided in this manual are for reference only, and the settings or data in the illustrations may not be exactly the same as the actual display you see on the product.

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Distributed by MyCardio LLC, 3200 E Cherry Creek South Dr, Suite 540, Denver, CO 80209, USA

Contact SleepImage: support@sleepimage.com

SleepImage Website: <https://sleepimage.com>

Trade Name: SleepImage Ring (Large)

Product Name: Pulse Oximeter

Model: PO2

Shenzhen Viatom Technology Co., Ltd.
901, Building West, Lepu Tower, No.66 Xingke Road, Xili Community, Xili Street, Nanshan District, Shenzhen, 518055, Guangdong, P.R. China
www.viatomtech.com

MedNet EC-REP GmbH
Borkstrasse 10, 48163 Muenster, Germany
Tel: +49 251 32266-0
Fax: +49 251 32266-22
Email: contact@mednet-ecrep.com

MediMap Ltd
2 The Drift, Thurston, Suffolk IP31 3RT, United Kingdom
Tel: +49 251 32266-0
Fax: +49 251 32266-22
Email: contact@mednet-ecrep.com

Australia Sponsor: SHARE INFO PTY LTD
4 Allnutt Court, Cheltenham, VIC 3192, Australia