SleepImage PO6 Fingertip

Model: PO6

User Manual

The SleepImage Fingertip is a version of the Viatom model PO6 and uses Bluetooth to communicate with the SleepImage Mobile App (App). The Fingertip must be paired with the SleepImage App to collect and transfer data to the SleepImage System, a Software as a Medical Device (SaMD).

Download the app by searching:

<https://sleepimage.com/patient-resources/>

The Fingertip displays Pulse Rate and SpO₂ data while in use. It should not be used to directly diagnose or treat any medical condition, and it does not store any data.

1. Introduction

1.1. Intended Use

This Pulse Oximeter is intended to be used for measuring, displaying, and storing adults' pulse oxygen saturation (SpO₂) and pulse rate in adults in home or healthcare facility environments. (EU AU)

The Pulse Oximeter is a non-invasive device intended for spot checking and/or continuous data collection of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (PR). This portable device is indicated for use in adult patients in home environments and clinical institutions except acute clinical environment. (US CA)

Intended Users: adults

1.2. Warnings and Cautions


- Do not use this device during MRI examination.
- Do not use this device with a defibrillator.
- Do not store the device in the following locations: direct sunlight; high temperatures or high levels of moisture, or heavy contamination; locations near to sources of water, fire, lint, or dust; or locations that are subject to strong electromagnetic influences.
- Do not use the device in a combustible environment.
- Never submerge the device in water or other liquids.
- Do not clean the device with acetone or other volatile solutions.
- Do not drop this device or subject it to strong impact.
- The device and accessories are provided non-sterile.
- Do not place this device in pressure vessels or gas sterilization device.
- Do not dismantle the device.
- Consult your doctor immediately if you experience symptoms that could indicate acute disease.
- Do not self-diagnose or self-medicate on the basis of this device without consulting your doctor. In particular, do not start taking any new medication or change the type and/or dosage of any existing medication without prior approval.
- Use only cables, sensors and other accessories specified in this manual.
- Prolonged continuous monitoring may increase the risk of undesirable changes in skin characteristics, such as irritation, reddening, blistering or burns.
- Do not open the device cover.

- The biocompatibility testing has been performed on the materials in contact with the person in accordance with ISO 10993.
- Do not place the device on a finger with edema or fragile tissue.
- Check the display values on the device before use. Do not use a device with no or unlikely values displayed.
- Check the sensor application site before going to sleep, during awakenings and at the end of the sleep test to determine correct positioning of the sensor and the circulation and skin sensitivity of the patient. Patient sensitivity varies depending on medical status or skin condition. This sensor may not be appropriate for patients with poor peripheral blood circulation or sensitive skin.
- The device has no alarm system.
- Local laws and regulations should be followed when disposing of the device and accessories.
- Do not wear the device while it is charging.
- Keep the charging cable away from children. It can cause strangulation.
- Keep the device out of reach of children and pets.
- This device is calibrated to display FUNCTIONAL OXYGEN SATURATION.
- Nail polish may affect the measuring accuracy, and too long fingernail may cause failure of measurement or inaccurate results.

1.3. Included in Box

- Device, User Manual, Warranty Card, Charging Cable.

2. Display Overview

3. Using the Device and App

3.1. Charging

Charge the battery before using by connecting the device to computer USB or USB charging adapter using the included USB cable. When fully charged, the device will power off automatically.

3.2. SleepImage App

The device uses Bluetooth to communicate with the SleepImage Mobile App. The device must be paired with the App to collect and transfer data to SleepImage.

Note: *DO NOT PAIR* in settings of your phone.

Follow the Mobile App Instructions for Use on www.sleepimage.com to download the App and pair the device *before* attempting data collection.

3.3. Power ON/OFF

Power ON:

The device will turn on automatically when worn.

Power OFF:

The device will turn off automatically after removal.

3.4. Start recording

Wear the device on index finger of left hand. The device will turn on automatically. **Start the recording in the App to begin data collection.**

Device fit and other notices:

- Wear the device on index finger of left hand.
- A proper fit is essential. Try to move the device to find the best fit.
- Avoid wearing the device loose on the finger, which causes inaccurate measure.
- Avoid excessive motion during recording.
- Avoid strong ambient light during recording.
- If the device is removed briefly during a recording session, simply put the device back on to resume recording.

3.6. Stop recording

After waking, Stop the recording in the App and remove the device which will automatically turn itself off.

3.7. Data Upload

After the recording is stopped in the App, the recording will upload to SleepImage. If study does not upload, check Troubleshooting in the SleepImage App Instructions for Use on www.sleepimage.com.

3.8 Cybersecurity

- The device can transmit clinical data (SpO₂ data) from the subject device to a wireless receiver for storage and display. The common receiver is APP software. So it is needed to install international mainstream anti-virus software such as AVG, Avast, McAfee, etc. In the running environment, even update the virus upgrade package in time and scan for viruses. For IOS system, the APP can run without antivirus software and firewalls because of the closed system and access control of the system. The end user needs to select the device ID before they want to conduct Bluetooth connection.
- In order to support the cybersecurity of this device, we provide information about the risk management with regards to cybersecurity which includes the assets, vulnerabilities, threats and controls. If you have any problems, please contact us.

Identification	Risk Description	Impact
Threats	CR1: Malicious attack such as Counterfeit command	Maybe the device failed to use and the treatment is delayed.
	CR2: Malicious attack such as Violent login attempt / password guessing attack	Maybe the treatment is delayed.
	CR3: Sensitive information including health data and the user information leakage	Maybe Tampering and Embezzlement of the user' information.
	CR4: Wrong operation of critical and/or sensitive information by authorized users	Maybe Tampering and Embezzlement of the user' information.

	CR5: Improper Authority Setting	Maybe the device failed to use.
Vulnerabilitie s	CR6: Access without authentication	Maybe the device failed to use.
	CR7: Sensitive information with encryption	Maybe the user' data is lost or stolen.
	CR8: Sensitive information with integrity	Maybe the user' data is lost or stolen.
Assets	CR9: Information Assets: Patient data; Device settings; Health data; Connection ID	Maybe the health data is revised or lost.
	CR10: Physical Assets :User interface; Network interface; Device assets (Operating system, Software etc.)	Maybe the device failed to use.

3) The device only includes the embedded software. The firmware is directly coded into the integrated system with on circuit unchangeable package. There is no software transferring during production or distribution. The manufacturer does the burning and upgrading of the subject device software and there is no case of being attacked by malware. Then installing antivirus software and firewalls in the device is unnecessary.

4) SleepImage regularly checks whether the released app software is the latest version, whether it can be downloaded and installed, and whether the software is attacked (MD5 hash value verification) to ensure software integrity and authenticity.

5) All updates and patches will be released to customer through distributor. A release note or letter will be issued and the version number will be changed accordingly.

6) Device software upgrade can be done remotely through our APP. When the device is connected to the APP, the APP will check the software version of the device, and when a new available version is found, an upgrade button will pop up and the user will authorise the upgrade operation.

7) Before discarding the device, it needs to be ensured that all health measurements and user account data stored in the device and the app are secure and confidential (e.g. deleted).

3.9. How to Check Battery level

Touch the power button to switch the display between showing readings, battery level and last four digits of the device's serial number (SN).

3.10. Unavailable Symbol

When this symbol displays on device screen, it indicates the readings is unavailable, which may be caused by:

- Excessive movement;
- Poor signal;
- Finger is too cold.

4. Maintenance

4.1. Time & Date

Device time automatically syncs from your mobile device.

4.2. Cleaning

Use a soft cloth moistened with water or alcohol to clean the device surface.

4.3. Recycling

Do not casually dispose the main unit and its accessories at the end of their service life. They shall be disposed in accordance with local regulations. When they are intended to be discarded, they must be sent to the appropriate facility for recycling.

5. Troubleshooting

Problem	Possible Cause	Possible Solution
Device does not turn on or no response	Battery may be low.	Charge battery and try again.
	Device may be damaged.	Please contact your local distributor.
The app cannot find or pair with the device	The device is off.	Turn on device.
	Bluetooth cannot work without location permission (Android)	Allow location access.
	The Bluetooth of your phone is off.	Turn on the Bluetooth in the phone.
	Your mobile device may not meet the minimum specifications required. (See www.sleepimage.com for information).	

6. Guide to Symbols

Symbol	Description
	Manufacturer
	Date of manufacture
	Serial number
	Indicates that the product should not be discarded as unsorted waste but must be sent to separate collection facilities for recovery and recycling.
	Follow Instructions for Use.
	Type BF Applied Part
	No alarm system
	MRI unsafe. Presents hazards in all MR environments as device contains strongly ferromagnetic materials.
IP22	Indicates that the product is protected against solid foreign objects of 12,5 mm Ø and greater; and protected against vertically falling water drops when enclosure tilted up to 15°.
	CE marking
	Medical device
	Authorized representative in the European community.
	UKCA marking.

	Authorized Representative in the United Kingdom.
	This product complies with the rules and regulations of the Federal Communication Commission.
	Prescription use
	Non-ionizing radiation.
	This device contains license-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada’s license-exempt RSS(s).
	Temperature limitation.
	Humidity limitation.
	Atmospheric pressure limitation.
	This product complies with verpackG.
	Our products and packaging can be recycled, don’t throw them away! Find where to drop them off on the www.quefairedemesdechets.fr site (Only applicable for French market)

7. Specifications

Environmental	Operating	Storage
Temperature	5 to 40°C	-10 to 45°C
Relative humidity (noncondensing)	10% to 80%	10% to 93%
Barometric	700 to 1060hPa	700 to 1060hPa
Protection against electric shock	Internally powered equipment	
Degree protection against electrical shock	Type BF	
Electro-magnetic compatibility	Group I, Class B	
Degree of dust & water resistance	IP22	
Weight	28 g	
Size	53mm×35mm×39 mm	
Battery	3.7Vd.c., Rechargeable Lithium-polymer	
Charge requirement	DC 5V, 1A	
Charge time	2-3 hours	
Battery life	12-14 hours for typical use	
Wireless	Bluetooth 4.2 BLE	
SpO ₂ range	0% to100%	
SpO ₂ accuracy (Arms)	70%-100%: ± 2%(Arms:1.88); 70%-80%: ± 3%; 80%-90%: ± 2%; 90%-100%: ± 2%; 0%-69%: not defined	
Pulse Rate range	30 to 250 bpm	

Pulse Rate accuracy	±2 bpm or ±2%, whichever is greater
A functional tester or SpO2 simulator can be used to determine the pulse rate accuracy.	
Wavelength / Max emission power	660nm/940nm, 0.8mW/1.2mW
Recorded parameters	SpO ₂ , Pulse Rate
Data storage	No data is stored on the device
Frequency range	2.402 – 2.480 GHz
Max RF power	-10 dBm
Expected service life	3 years

8. FCC Statement

FCC Warning:
FCC ID: 2AD XK-4708
Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and
(2) this device must accept any interference received, including interference that may cause undesired operation.
Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
-Reorient or relocate the receiving antenna.
-Increase the separation between the equipment and receiver.
-Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
-Consult the dealer or an experienced radio/TV technician for help.
The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

9. IC Caution

IC:29845-S001
This device contains license-exempt transmitter(s)/ receiver(s) that comply with Innovation, Science and Economic Development Canada’s license-exempt RSS(s). Operation is subject to the following two conditions:
(1) This device may not cause interference.
(2) This device must accept any interference, including interference that may cause undesired operation of the device.
Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

10. Electromagnetic Compatibility

The device meets the requirements of IEC 60601-1-2: 2020.

Shenzhen Viatom Technology Co., Ltd.
901, Building West, Lepu Tower, No.66 Xingke Road, Xili Community, Xili Street, Nanshan District, Shenzhen, 518055, Guangdong, P.R. China
www.viatomtech.com

MedNet EC-REP GmbH
Borkstrasse 10, 48163 Muenster, Germany

MediMap Ltd
2 The Drift, Thurston, Suffolk IP31 3RT, United Kingdom

Australia Sponsor: SHARE INFO PTY LTD
Add: 4 Allnutt Court, Cheltenham, VIC 3192, Australia

Trade name: SleepImage PO6 Fingertip
Product name: Pulse Oximeter Model: PO6
Version: V1.1 Date:August 21, 2026 PN:255-04060-C3

Distributed by MyCardio LLC, 3200 E Cherry Creek South Drive, Suite 540, Denver, CO 80209, USA

SleepImage Website: www.sleepimage.com
Contact SleepImage: support@sleepimage.com

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